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2. MSH - Message header segment

The MSH segment defines the intent, source, destination, and some specifics of the syntax of a message

MSH field definitions

SEQ	HL7 LEN	HL7 DT	HL7 R/O	ELEMENT NAME	ELR Usage
1	4	ST	R	Encoding characters	^~\&
2	180	HD	O	Machine Name	MachineName
3	180	HD	O	Company Name	CompanyName
4	180	HD	O	Receiving application	HIS
5	180	HD	O	Receiving facility	BMS-HOSxP
6	26	TS	R	Date/Time of message	Create Message Date Time (yyyymmddhhnnss)
7	40	ST	O	Security	
8	7	CM	R	Message type/Trigger Event	ORM^O01, ORU^R01, ADT^A40, ADT^08
9	20	ST	R	Message control ID	Generate file by file
10	3	PT	R	Processing ID	Set to 'P'
11	60	ID	R	Version ID	2.3

3. PID - Patient Identification

The PID segment is used by all applications as the primary means of communicating patient identification information. This segment contains permanent patient identifying and demographic information that, for the most part, is not likely to change frequently. Usage notes: PID patient identification The assigning facility ID, the fourth component of the patient identifiers, is a string of up to six characters that is uniquely associated with the facility that originally assigned the number.

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4. PV1 - Patient visit

The PV1 segment is used by Registration/Patient Administration applications to communicate information on a visit-specific basis. This segment can be used to send multiple-visit statistic records to the same patient account or single-visit records to more than one account.

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5. OBR - Observation request segment

The Observation Request (OBR) segment is used to transmit information specific to an order for a diagnostic study or observation, physical exam, or assessment. The Observation Request segment defines the attributes of a particular request for diagnostic services (e.g., laboratory, EKG) or clinical observations (e.g., vital signs or physical exam).

OBR field definitions

SEQ	HL7 LEN	HL7 DT	HL7 R/O	ELEMENT NAME	ELR Usage
1	4	SI	C	Set ID - Observation Request	1 (id of data)
2	75	EI	C	Placer Order Number	**Not Used
3	75	EI	C	Filler Order Number	**Not Used
4	200	CE	R	Universal Service Identifier	**Not Used
5	2	ID	B	Priority	**Not Used
6	26	TS	B	Requested Date/Time	20200429075741
7	26	TS	C	Observation Date/Time	**Not Used
8	26	TS	O	Observation End Date/Time	**Not Used
9	20	CQ	O	Collection Volume	**Not Used
10	60	XCN	O	Collector Identifier	**Not Used
11	1	ID	O	Specimen Action Code	**Not Used
12	60	CE	O	Danger Code	**Not Used
13	300	ST	O	Relevant Clinical Information	**Not Used
14	26	TS	C	Specimen Received Date/Time	20200429075741
15	300	CM_SPS	O	Specimen Source	**Not Used

6. Observation result Segment (OBX)

The OBX segment is used to transmit a single observation or observation fragment. It represents the smallest indivisible unit of a report. Its principal mission is to carry information about observations in report messages. But the OBX can also be part of an observation order. In this case, the OBX carries clinical information needed by the filler to interpret the observation the filler makes. OBX is also found in other HL7 messages that need to include patient clinical information.

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